


Dil/DiA axonal tracing

1. Brains are dissected in PBS, and post-fixed with immersion in 4% PFA at least overnight. (Labeling can be done after one week in PFA)
2. Determine the dye placement location under a dissecting scope. You can first take a digital image of the brain and superimpose a pre-constructed grid on the image to determine the location.
3. Make a small hole in the pia at the insertion site using a 36-gauge syringe needle.
4. Use insert pins (No. 26007-02 and No 26018-17 or 26016-12, FST) to insert a single Dil or DiA crystal (D3911 and D3883, Invitrogen) into the brain tissue. Make sure to use the same pine for a particular dye. It is also important not to penetrate the crystal into the underlying white matter.
5. Image the brain to record dye placement.
6. Keep brain in 4% PFA at 37°C. The time of incubation depends on the stage of the brain (see following table).
7. Once diffusion is complete, the brains are sectioned with a vibratome at 100 μm thickness.
8. Sections are mounted with 30% sucrose. (Don't use commercial mounting media, which often contain glycerol). Make sure no detergent is present during sectioning and mounting.
9. Take picture as soon as possible.

Table I (based on Gurung & Fritsch, *J Comp Neurol*, 2004)

Age	E13.5	E14.5	E16.5	E18.5	P0	>P2
Days	8	10	10-13	14-15	18	4-6 weeks



A. Schematic representation of topographic connections between specific cortical areas with thalamic nuclei. B. Insertions of DiA and Dil crystals. C. Retrograde labeling of thalamic neurons.