

Preparation of electrocompetent cells

(Based on Protocol 26 from Molecular Cloning)

Day One: This step can be done at the end of a day.

1. Inoculate a single colony of *E. coli* from a fresh agar plate into a flask containing 50 ml of LB medium. Incubate the culture overnight at 37°C in a rotary shaker (250 rpm).

Day two: It starts in the morning and takes about ~ hours. Have these solutions ready: prewarmed 1 liter LB medium, pre-cooled 500 ml of Ultra-pure water, 10% glycerol, and GYT medium.

2. Inoculate two aliquots of 500 ml of prewarmed LB medium in separate 2-liter flasks with 25 ml of the overnight bacterial culture. Incubate the flasks at 37°C with agitation (300 rpm in a rotary shaker). Measure the OD₆₀₀ of the growing bacterial cultures every 20 min.
This density is usually achieved after ~2.5 hours of incubation for DH5a (or DH10B?). It is essential that the density does not exceed OD₆₀₀, 0.4.
3. When the OD₆₀₀ of the cultures reaches 0.4, rapidly transfer the flasks to an ice-water bath for 15-30 minutes. Swirl the culture occasionally to ensure that cooling occurs evenly. In preparation of the next step, place the centrifuge bottles in an ice-water bath, and pre-cool the rotor and chamber of centrifuge to 4°C.
For maximum efficiency of transformation, it is crucial that the temperature of the bacteria not rise above 4°C at any stage in the protocol.
4. Transfer the cultures to ice-cold centrifuge bottles (4X 250 ml). Harvest the cells by centrifugation at 1,000g (2,500 rpm in a Sorvall SLA 1500 rotor) for 15 min at 4°C. Decant the supernatant, resuspend the cell pellet in 500 ml of ice-cold pure H₂O (2X 250 ml bottle).
5. Harvest the cells by centrifugation at 1,000g for 20 min at 4°C. Decant the supernatant and resuspend the cell pellet in 250 ml of ice-cold 10% glycerol (1X250 ml bottle).
This should be done immediate after the centrifugation finishes. Take care when decanting because the bacterial pellets lose adherence in 10% glycerol.
6. Harvest the cells by centrifugation at 1,000g for 20 min at 4°C. Decant the supernatant and resuspend the cell pellet in 10 ml of ice-cold 10% glycerol. (Cells can be transferred to 50 ml Falcon tube and spun in a bench-top centrifuge).
7. Harvest the cells by centrifugation at 1,000g for 20 min at 4°C. Carefully decant the supernatant and use a Pasteur pipette attached to a vacuum line to remove any remaining drops of buffer. Resuspend the pellet in 1 ml of ice-cold GYT medium.
This is best done by gentle swirling rather than pipetting or vortexing.
8. Measure the OD₆₀₀ of a 1:100 dilution of the cell suspension. Dilute the cell suspension to a concentration of 2 X 10¹⁰ to 3 x 10¹⁰ cells/ml (1.0 OD₆₀₀=~2.5 X 10⁸ cells/ml) with ice-cold GYT medium.

9. Transfer 40 μ l of the suspension to an ice-chilled electroporation cuvette (0.2-cm gap) and test whether arching occurs when an electrical discharge is applied. If so, wash the remainder of the cell suspension once more with ice-cold GYT medium to ensure that the conductivity of the bacterial suspension is sufficiently low.
10. Dispense 40- μ l aliquots of the cell suspension into sterile, ice-cold 0.5-ml microfuge tubes, drop into a bath of liquid nitrogen, and transfer to a -70°C freezer.

Solutions:

Glycerol (10% v/v)

GYT Medium

10% (v/v) glycerol

0.12% (w/v) yeast extract

0.25% (w/v) tryptone

Sterilize the medium by passing it through a prerinsed 0.22- μm filter. Store in 2.5 ml aliquots at 4°C .