

## Nuclear and Cytoplasmic Protein Extraction (Thermo Scientific)

### Cell Culture Preparation

1. For adherent cells, harvest with trypsin-EDTA and then centrifuge at  $500 \times g$  (250 rpm) for 5 minutes. For suspension cells, harvest by centrifuging at  $500 \times g$  for 5 minutes.
2. Wash cells by suspending the cell pellet with PBS.
3. Transfer  $1-10 \times 10^6$  cells to a 1.5 ml microcentrifuge tube and pellet by centrifugation at  $500 \times g$  for 2-3 minutes.
4. Use a pipette to carefully remove and discard the supernatant, leaving the cell pellet as dry as possible.
5. Add ice-cold CER I (add protease inhibitors 1x before use) to the cell pellet (Table 1). Proceed to Cytoplasmic and Nuclear Protein Extraction, using the reagent volumes indicated in Table 1.

Packed Cell Volume ( $\mu$ l)	CER I ( $\mu$ l)	CER II ( $\mu$ l)	NER ( $\mu$ l)
10	100	5.5	50
20	200	11	100
50	500	27.5	250
100	1,000	55	500

\*For HeLa cells,  $2 \times 10^6$  cells is equivalent to 20  $\mu$ l packed cell volume.

### Cytoplasmic and Nuclear Protein Extraction

**Note:** Scale this protocol depending on the cell pellet volume (Tables 1 and 2). Maintain the volume ratio of CER I:CER II:NER reagents at 200:11:100  $\mu$ l, respectively.

1. Vortex the tube vigorously on the highest setting for 15 seconds to fully suspend the cell pellet. Incubate the tube on ice for 10 minutes.
2. Add ice-cold CER II to the tube.
3. Vortex the tube for 5 seconds on the highest setting. Incubate tube on ice for 1 minute.
4. Vortex the tube for 5 seconds on the highest setting. Centrifuge the tube for 5 minutes at maximum speed in a microcentrifuge ( $\sim 16,000 \times g$ ).
5. Immediately transfer the supernatant (cytoplasmic extract) to a clean pre-chilled tube. Place this tube on ice until use or storage (see Step 10).
6. Suspend the insoluble (pellet) fraction produced in Step 4, which contains nuclei, in ice-cold NER (add protease inhibitors 1x before use).
7. Vortex on the highest setting for 15 seconds. Place the sample on ice and continue vortexing for 15 seconds every 10 minutes, for a total of 40 minutes.
8. Centrifuge the tube at maximum speed ( $\sim 16,000 \times g$ ) in a microcentrifuge for 10 minutes.
9. Immediately transfer the supernatant (nuclear extract) fraction to a clean pre-chilled tube. Place on ice.
10. Store extracts at  $-80^\circ\text{C}$  until use.