

PCR genotyping

Preparation of DNA lysates

1. Cut tail (3~5 mm).
2. Prepare lysis buffer. To each sample add 50 μ l of lysis buffer with fresh ProK to a final concentration of 200 μ g/ml (1.5 μ l ProK each). Incubate at 60°C overnight.
3. Inactivate ProK activity by heating lysates at 100°C for 10 min.
4. Quick spin for 10"

Amount of lysis buffer:	<u>Yolk sac</u>
Tail or toe: 100 μ l	E8.5: 15 μ l
E7.5 embryos: 15 μ l	E9.5-E10.5: 50 μ l
E8.5 embryo: 30 μ l	>E10.5: 100 μ l

PCR reaction

1. Prepare PCR reaction mixture and aliquot 19 μ l to each tube.
2. Add 1 μ l of DNA lysate as template.
3. Put tubes in PCR machine and initiate designed program.

Note. To have a more robust and cleaner PCR amplification, I routinely use "hot start" PCR. To do this, prepare the reaction on ice, start PCR program and wait until the block reach 94°C to put in tubes quickly.

Run gel to analyze PCR amplified products.

Lysis buffer I

50 mM Tris pH8.8
1mM EDTA (pH8.0)
0.5% Tween 20
Add protease K 200 μ g/ml to the final concentration

Lysis buffer II

50 mM KCl	1M	5 ml
10 mM Tris pH8.3	1M	1 ml
2.5 mM MgCl ₂	1M	0.25 ml
0.1 mg/ml Gelatin	1 mg/ml	10 ml
0.45% NP40	100%	0.45 ml
0.45% Tween 20	100%	0.45 ml

Add proteinase K to 100 μ g/ml. 400 μ l for tail, 100 μ l E8.5 YS, 200 μ l E9.5 YS.

PCR reaction

Gbx2/Neo PCR

	1X	8X	16X	24X	32X
Buffer	2	16	32	48	64
H ₂ O	15.2	121.6	243.2	364.8	486.4
Gbx2-R (20 µM)	0.5	4	8	12	16
Gbx2-F (20 µM)	0.5	4	8	12	16
Neo (15-) (1µg/µl)	0.2	1.6	3.2	4.8	6.4
Neo (255-) (1µg/µl)	0.2	1.6	3.2	4.8	6.4
dNTPs (25 mM)	0.2	1.6	3.2	4.8	6.4
Taq	0.2	1.6	3.2	4.8	6.4
Total	19	152	304	456	608

Neo: 250 bp
Gbx2: 180 bp
Gbx2-flox: 220 bp
Run 2~3% gel

Neo (15-36): GAACAAGATGGATTGCACGCAG
Neo (266-245): TTCAGTGACAACGTCGAGCACA
Gbx2 (1586-1565): TGCTTGGATGTCCACATCTAGG
Gbx2 (1404-1425): CTGTTACGTTAGCAGGTTTCGC

LacZ/RAP PCR

	1X	8X	16X	24X	32X
Buffer	2	16	32	48	64
H ₂ O	15	120	240	360	480
LacZ1 (20 µM)	0.5	4	8	12	16
LacZ2 (20 µM)	0.5	4	8	12	16
RAP1 (20 µM)	0.3	2.4	4.8	7.2	9.6
RAP2 (20 µM)	0.3	2.4	4.8	7.2	9.6
dNTPs (25 mM)	0.2	1.6	3.2	4.8	6.4
Taq	0.2	1.6	3.2	4.8	6.4
Total	19	152	304	456	608

LacZ: 800 bp
RAP: 600 bp
Run 1% gel

Program /JAMES/GBX2

1. 95°C 3'
2. 95°C 55"
3. 58°C 50"
4. 72°C 1'
5. Go to step 2 and repeat 30 times
6. 4°C for ever
7. End

Program /JAMES/LACZ

1. 95°C 3'
2. 95°C 40"
3. 72°C 45"
4. Go to step 2 and repeat 30 times
5. 4°C for ever
6. End

Otx1/Otx2 PCR

	1X	8X	16X	24X	32X
Buffer	2	16	32	48	64
H ₂ O	16.2	129.6	259.2	388.8	518.4
Otx2-r (0.5µg/µl)	0.1	0.8	1.6	2.4	3.2
Otx2-f (0.5µg/µl)	0.1	0.8	1.6	2.4	3.2
Otx1-r (1µg/µl)	0.1	0.8	1.6	2.4	3.2
Otx1-f (1µg/µl)	0.1	0.8	1.6	2.4	3.2
dNTPs (25 mM)	0.2	1.6	3.2	4.8	6.4
Taq	0.2	1.6	3.2	4.8	6.4
Total	19	152	304	456	608

Otx2: 220 bp
hOtx1: 320 bp
Run 2~3% gel

Program /JAMES/OTX

1. 95°C 3'
2. 95°C 50"
3. 64°C 55"
4. 72°C 1'
5. Go to step 2 and repeat 30 times
6. 4°C for ever
7. End

Gbx2/Neo-Promoter

	1X	8X	16X	24X	32X
Buffer	2	16	32	48	64
H ₂ O	15.4	124	247	370	493
Gbx2-R (20 µM)	0.5	4	8	12	16
Gbx2-F (20 µM)	0.5	4	8	12	16
Neo-PrmR (1µg/µl)	0.2	0.8	1.6	2.4	3.2
dNTPs (25 mM)	0.2	1.6	3.2	4.8	6.4
Taq	0.2	1.6	3.2	4.8	6.4
Total	19	152	304	456	608

Gbx2: 180 bp
Neo: 300 bp
Run 2% gel

Program /JAMES/Gbx2

1. 95°C 3'
2. 95°C 50"
3. 59°C 55"
4. 72°C 1'
5. Go to step 2 and repeat 30 times
6. 4°C for ever
7. End

Flp PCR

	1X	8X	16X	24X	32X
Buffer	2	16	32	48	64
H ₂ O	16.2	129.6	259.2	388.8	518.4
SD41 (20μM)	0.2	1.6	3.2	4.8	6.4
SD42 (20 μM)	0.2	1.6	3.2	4.8	6.4
dNTPs (25 mM)	0.2	1.6	3.2	4.8	6.4
Taq	0.2	1.6	3.2	4.8	6.4
Total	19	152	304	456	608

Flp: 750 bp

Run 1% gel

SD41: GTGGATCGATCCTACCCCTTGCG

SD42: GGTCCAAGTGCAGCCCAAGCTTCC

Cre PCR

	1X	8X	16X	24X	32X
Buffer	2	16	32	48	64
H ₂ O	16.2	129.6	259.2	388.8	518.4
Cre-f (20 μM)	0.2	1.6	3.2	4.8	6.4
Cre-r (20 μM)	0.2	1.6	3.2	4.8	6.4
dNTPs (25 mM)	0.2	1.6	3.2	4.8	6.4
Taq	0.2	1.6	3.2	4.8	6.4
Total	19	152	304	456	608

Cre: ~300 bp

Run 1~2% gel

CreF: TAAAGATATCTCACGTAAGTACGCGTG

CreR: TCTCTGACCAGAGTCATCCTTAGC

Gbx2-cDNA PCR

	1X	8X	16X	24X	32X
Buffer	2	16	32	48	64
H ₂ O	16.2	129.6	259.2	388.8	518.4
Gbx2 B2 (20 μM)	0.2	1.6	3.2	4.8	6.4
Gbx2-5a' (20 μM)	0.2	1.6	3.2	4.8	6.4
dNTPs (25 mM)	0.2	1.6	3.2	4.8	6.4
Taq	0.2	1.6	3.2	4.8	6.4
Total	19	152	304	456	608

Gbx2 B2: CAGTGGAAATCTTTCTCC

Gbx2-5a': AGACGGCAAAGCCTTCTTGG

Program /JAMES/FLP

1. 95°C 3'
2. 94°C 55"
3. 67°C 55"
4. 72°C 1'
5. Go to step 2 and repeat 30 times
6. 4°C for ever

Program /JAMES/CRE

1. 95°C 2'
2. 95°C 40"
3. 59°C 1'
4. 72°C 50"
5. Go to step 2 and repeat 29 times
6. 4°C for ever
7. End

Gbx2 cDNA: ~400 bp

Run 1~2% gel

Program /JAMES/GBX2

1. 95°C 3'
2. 94°C 45"
3. 55°C 50"
4. 72°C 50"
5. Go to step 2 and repeat 25 times
6. 4°C for ever
7. End

Gbx2/Cre PCR

	1X	8X	16X	24X	32X
Buffer	2	16	32	48	64
H ₂ O	15.2	121.6	243.2	364.8	486.4
Gbx2-R (20 µM)	0.5	4	8	12	16
Gbx2-F (20 µM)	0.5	4	8	12	16
Cre-f (1 µg/µl)	0.2	1.6	3.2	4.8	6.4
Cre-r (1 µg/µl)	0.2	1.6	3.2	4.8	6.4
dNTPs (25 mM)	0.2	1.6	3.2	4.8	6.4
Taq	0.2	1.6	3.2	4.8	6.4
Total	19	152	304	456	608

Neo: 250 bp
Gbx2: 180 bp
Gbx2-flox: 220 bp
Cre: 700 bp
Run 2~3% gel

Program /JAMES/CRE

1. 95°C 2'
2. 95°C 40"
3. 59°C 1'
4. 72°C 50"
5. Go to step 2 and repeat 29 times
6. 4°C for ever
7. End

Gbx2-Neo-promoter/Cre

	1X	8X	16X	24X	32X
Buffer	2	16	32	48	64
H ₂ O	15.5	124	248	372	496
Gbx2-R (20 µM)	0.5	4	8	12	16
Neo-PrmR (1 µg/µl)	0.2	1.6	3.2	4.8	6.4
Cre-f (1 µg/µl)	0.2	1.6	3.2	4.8	6.4
Cre-r (1 µg/µl)	0.2	1.6	3.2	4.8	6.4
dNTPs (25 mM)	0.2	1.6	3.2	4.8	6.4
Taq	0.2	1.6	3.2	4.8	6.4
Total	19	152	304	456	608

LacZ: 800 bp
RAP: 600 bp
Run 1% gel

Program /JAMES/CRE

1. 95°C 2'
2. 95°C 40"
3. 59°C 1'
4. 72°C 50"
5. Go to step 2 and repeat 29 times
6. 4°C for ever
7. End

ROSA26 locus

	1X	8X	16X	24X	32X
Buffer	2	16	32	48	64
H ₂ O	17.3	138.2	276.8	415.2	553.6
RR1/2/3 (1µg/µl)	0.3	2.4	4.8	7.2	9.6
dNTPs (25 mM)	0.2	1.6	3.2	4.8	6.4
Taq	0.2	1.6	3.2	4.8	6.4
Total	19	152	304	456	608

Program /JAMES/ROSA.

1. 94°C 2'
2. 95°C 1'
3. 61°C 1'
4. 72°C 1'
5. Go to step 2 and repeat 29 times
6. 72°C 10'
7. 4°C for ever
8. End

Wt band: 580 bp

Mutant: 254 bp

Run 1% gel

RR1: AAAGTCGCTCTGAGTTGTTAT
RR2: GCGAAGAGTTTGTCTCAACC
RR3: GGAGCGGGAGAAATGGATATG

