

Ten-minute plasmid DNA minipreps

(Based on Zhou, C., Yujun, Y and Jong, A.Y. (1990) BioTechniques 8; 172-173)

1. Centrifuge ~1.2 ml culture^a 20 sec in a 1.5 ml tube.
2. Remove supernatant by aspiration. Add 50 μ l of Solution I^b; vortex at highest speed to resuspend pellet.
3. Add 300 μ l TENS^c, vortex 2-5 seconds (it becomes sticky)^d.
4. Add 150 μ l 3M NaOAc, pH5.2; invert 10~16 times.
5. Centrifuge at top speed for 5 minutes^d. At the meantime, add 900 μ l of 100% EtOH to new tubes.
6. Pour supernatant to the tube with EtOH, avoiding the precipitate. Mix by vortexing.
7. Centrifuge at top speed for 5 minutes^f.
8. Pour off supernatant. Add 600 μ l of 70% EtOH (you can dispense 70% EtOH from a bottle). Pour off 70% EtOH, and quick spin
9. Remove most 70% EtOH; Air dry pellet briefly.
10. Resuspend in 20-40 μ l TE^g. This usually results in ~ 100 ng/ μ l.

Notes:

- a. Using sterile technique, pick a colony into 1.5 ml LB broth in a 15 ml Falcon (2059) tube. Shake vigorously at 37°C overnight or a few hours until cloudy.
- b. Solution I (50 mM glucose, 25 mM Tris, 10 mM EDTA, pH8.0. For 200 ml solution, 10 ml of 1M glucose (or 1.8 g powder), 5 ml of 1M Tris.HCl, pH7.5, 4 ml of 0.5M EDTA, pH8.0, add H₂O to final 200 ml, add RNaseA to final concentration 100 μ g/ml)
- c. TENS=TE with 0.1N NaOH and 0.5% SDS (0.5 ml 10N NaOH, 1.25 ml 20% SDS for 50 ml solution). Store at RT.
- d. Keep on ice to prevent chromosomeal degradation if you need to take more than 10 minutes before step 4.
- e. This pellets cell debris and chromosomal DNA.
- f. This pellets plasmid DNA and RNA.
- g. Plasmid yield is 2.3 μ g from a 1.5 ml culture. Usually we add 30 μ l. Add more or less depending upon the size of the pellet. Use 1 μ l for digestion.