

## Retrieving DNA by Gap Repairing

Day One: This step can be done at the end of a day.

1. Inoculate a single colony of EL350 cells containing a BAC from a fresh agar plate into a 14-ml tube containing 5 ml of LB medium. Incubate the culture overnight at **32°C** in a rotary shaker (250 rpm).

Day two: It starts in the morning and takes about ~ hours. Have these solutions ready: prewarmed 1 liter LB medium, pre-cooled 500 ml of Ultra-pure water, 10% glycerol, and GYT medium.

2. Inoculate 20 ml of prewarmed LB medium in a 500-ml flask with 1 ml of the overnight bacterial culture ( $OD_{600} \sim 1.2$ ). Incubate the flasks at **32°C** with agitation (180 rpm in a rotary shaker). Measure the  $OD_{600}$  of the growing bacterial cultures every 20 min.

*This density is usually achieved after ~2.5 hours of incubation for DH5a (or DH10B?). It is essential that the density does not exceed  $OD_{600}$  0.4.*

3. When the  $OD_{600}$  of the cultures reaches 0.5 (~2 hr), incubate the flask at **42°C** for 15 min.
4. Immediately transfer the flasks to an ice-water bath, swirl the culture occasionally to ensure that the temperature drops as fast as possible. Incubate on ice for an additional 5 min. In preparation of the next step, place the centrifuge bottles in an ice-water bath, and pre-cool the rotor and chamber of centrifuge to 4°C.  
*For maximum efficiency of transformation, it is crucial that the temperature of the bacteria not rise above 4°C at any stage in the protocol.*
5. Transfer the cultures to an ice-cold 50-ml Falcon tube. Harvest the cells by centrifugation at 1,000g for 20 min at 0°C in a bench-top centrifuge. Decant the supernatant, resuspend the cell pellet in 20 ml of ice-cold pure H<sub>2</sub>O.
6. Harvest the cells by centrifugation at 1,000g for 10 min at 0°C. Decant the supernatant and resuspend the cell pellet in 10 ml of ice-cold 10% glycerol.  
*This should be done immediate after the centrifugation finishes. Take care when decanting because the bacterial pellets lose adherence in 10% glycerol.*
7. Harvest the cells by centrifugation at 1,000g for 20 min at 0°C. Decant the supernatant and resuspend the cell pellet in 2 ml of ice-cold 10% glycerol, transfer to two microfuge tubes.
8. Harvest the cells by centrifugation at 1,000g for 5 min at 0°C. Carefully decant the supernatant and use a Pasteur pipette attached to a vacuum line to remove any remaining drops of buffer. Resuspend the pellet in 50  $\mu$ l of ice-cold GYT medium.  
*This is best done by gentle swirling rather than pipetting or vortexing.*
9. Add 1-2  $\mu$ l of linearized retrieving vector (200-500 ng) to the cells and electroporate.

Solutions:

Glycerol (10% v/v)

GYT Medium

10% (v/v) glycerol

0.12% (w/v) yeast extract

0.25% (w/v) tryptone

Sterilize the medium by passing it through a prerinsed 0.22- $\mu$ m filter. Store in 2.5 ml aliquots at 4°C.