

PRIORITIES

DAILY PRIORITY



COVID-19 Symptoms
 If you are experiencing any symptoms, please do not come to work unless you have been cleared by a physician.

Possible COVID-19 Symptoms

- Fever
- Dry Cough
- Body Aches
- Loss of smell or taste
- Headache



UConn Health COVID-19 Call Center
 860-679-3199

Cleaning Protocol



Please clean prior to starting your microscope session.

Do not start the equipment until you have thoroughly cleaned the work area and equipment. You do not want to be charged for your time spent cleaning.

What needs to be cleaned?

- If you are going to touch it, then clean it.
- Wear a pair of gloves during the cleaning session.
- It is highly recommended to continue to wear gloves during your session, but please change your gloves after the initial cleaning.



CLEANING AGENTS



Disinfectant Wipes or Kimwipes Saturated with 70% Ethanol



70% Ethanol in Spray Bottle

- Light Switches
- Component Switches
- Incubator Handles and Knobs
- Microscope Components
- Immersion oil bottle
- Work surfaces

Remember:
 Allow the cleaning solution to be in contact with the item for 3 minutes or allow it to dry.

- Light Switches
- Component Switches
- Incubator Handles and Knobs
- Vinyl Chair



Throw your waste away!



Soap and water works best; hand sanitizer when there isn't a sink around.



Don't have others clean up after you.

PLEASE

DO

- Allow for 30 minutes between bookings for cleaning time.
- Continue to wear a mask even when you are alone in the room.
- Ask for assistance in advance. Things can go wrong but if you need help, plan for it.

DON'T

- Have multiple people in the lab.
- Remove your mask because you are alone; it is a multi-user room.
- Use the microscope if you have not scheduled your time.