

Prenatal Screening Laboratory

Test name	Specimen type	Notes
First trimester maternal serum screen (PAPP-A & hCG)	SERUM: one lg. Red or one lg. SST	Red must be spun and serum drawn off. 10-13.9 weeks.
Second trimester maternal serum (Quad) screen	SERUM: One lg. Red or one lg. SST	Red must be spun and serum drawn off. 14-21.9 weeks.
AFP only.	SERUM: One lg. Red or one lg. SST	Red must be spun and serum drawn off. 15-23.9 weeks.